

Methods for freezing mouse sperm in urgent situation

~Easy & Quick (EQ) method~

RIKEN Bioresource Research Center
Integrative Developmental Engineering Division
Updated technologies #1

May 1st, 2025

Contents

- **Introduction**
- **Reagents and Equipment**
- **Methods**
 1. Preparation of freezing solution, microtubes/tube holders
 2. Collection of cauda epididymis and freezing sperm
 3. Preparation of preincubation and fertilization dishes for IVF
 4. Thawing of frozen sperm and IVF
- **Results**
- **References**
- **Appendix**
- **Contact**

Introduction

Advantages

- Effective method for preserving sperm from many mice in case of disaster or emergency, and can be used for IVF after thawing.
- No liquid nitrogen required and can be stored in the freezer at -80°C .
- Can be performed quickly, even by inexperienced users, and can preserve sperm from 100 males in a short time.

Key points

- Simple method: Place the epididymis in a microtube with a 20% raffinose solution, cut with scissors and store directly in the -80°C freezer.
- If long-term preservation and high fertilization rates are required, standard methods of cryopreservation are preferable.

Reagents and Equipment

Raffinose: Becton Dickinson, cat#217410

Microtube (1.5 mL)

Aluminum foil

Fine scissors


Fine forceps

Tube storage boxes.

Heating devices (e.g. water bath) to dissolve at 50°C


Methods

1. Preparation of freezing solution, microtubes/tube stands




Methods

2. Collection of cauda epididymis and sperm freezing



Methods


3. Preparation of preincubation and fertilization dishes for IVF



***For more information on culture media and IVF conditions, see reference**

Methods

4. Thawing of frozen sperm and performing of IVF



Results

1. Fertilization rates using sperm frozen by standard and EQ methods.

Strain of males	Method	Ca ²⁺ concentration (mM)	No. of epididymides used	Total no. (%) of oocytes fertilized/inseminated	Mean % of fertilization (\pm SEM)
B6N	Standard	2	8	427/469 (91.0)	91.4 (\pm 1.5)
	EQ	2	7	209/401 (52.1)	49.9 (\pm 7.0)**
		5	7	176/307 (57.3)	56.7 (\pm 4.9)**
B6J	Standard	2	10	583/754 (77.3)	77.3 (\pm 1.4)
	EQ	2	7	118/380 (31.1)	30.8 (\pm 2.0)**
		5	7	150/301 (49.3)	50.9 (\pm 10.5)*

Result:

Even with the simple EQ method,
fertilization rates of about 50% have been achieved.

Results

2. Full-term development of fertilized embryos with sperm stored for one month using the EQ method.

Strain of males	No. (%) of females used			No. (%) of embryos				
	Transferred	Pregnant	(%)	Transferred	Implanted	(%)	Developed to offspring	(%)
B6N	3	3	(100)	45	39	(87)	26	(58)
B6J	4	4	(100)	60	49	(82)	33	(55)

Result:

Offspring were obtained at high developmental rates (>50%) in both strains.

Reference

- Mochida K, Hasegawa A, Shikata D, Itami N, Hada M, Watanabe N, Tomishima T, Ogura A. [Easy and quick \(EQ\) sperm freezing method for urgent preservation of mouse strains](https://doi.org/10.1038/s41598-021-93604-y). Sci. Rep. 2021, 11:14149.
<https://doi.org/10.1038/s41598-021-93604-y>

Appendix

Tips for sperm freezing

- It is important to cool the tube quickly. First, an aluminium foil holder should be made to fit the shape of the tube.
- This simple method can be used to preserve the sperm of a large number of mice in a short time. It is better to store in LN₂ as soon as possible rather than in the freezer for a longer period to reduce freezing damage.
- Under usual situation, we recommend using standard freezing methods for IVF that are easy to use.
- Sperm cannot be preserved at temperatures beyond -80°C.

Contact

- Integrative Developmental Engineering Division

Keiji Mochida : keiji.mochida@riken.jp

Kohtaro Morita : kohtaro.morita@riken.jp